

RNA Isolation (TRIZOL method)

- 01. Wash sample with PBS or MQ.**
- 02. Homogenize sample in 1ml Trizol (50-100mg sample)**
- 03. R.T. incubation for 10min and 14K for 10min at 4 °C.**
- 04. Add 200 µl chloroform to supernatant (1/5 volume of trizol).**
- 05. R.T. incubation for 3-5min and 12k for 15min at 4 °C.**
- 06. Add 1 volume PCI to supernatant, 15k for 5min at 4 °C.**
- 07. Add 1 volume (500 µl) isopropanol to supernatant.**
- 08. R.T. 10-15min and 12k for 10-15min at 4 °C.**
- 09. Dry and redissolve in 400 µl DEPC MQ.**
- 10. Add 1ml (2.5volume) 100% EtOH and 40 µl 3M NaOAc.**
- 11. 14K for 10min at 4 °C.**
- 12. Wash with 75% EtOH.**
- 13. Air dry pellet for 10min. and dissolve in 20-50 µl DEPC water.**
- 14. Heat for 10min at 60 °C to help dissolve RNA.**

6, 9-11, 14 are optional.